

SinoAngio series

DSA Injectors

SinoAngio



Perfect display quality



- 12.1 inch color LCD control display, large visual angle, graphic interface

Specifications of DSA High-Pressure Injector

Model	SinoAngio-1200	SinoAngio-600
Syringe	Single syringe	
Size of syringe	150mL sterile syringe	
Injection volume	0.1mL ~ Volume of syringe	
Injection delay	0~600 s	
Injection phases	1~8 phases	
Injection protocol	Up to 120 protocols	
Syringe heater	35°C±5°C (≥40min)	
Tilt lockout	Rotate head to arm-down, minimize risk of air embolish	
Control console	12.1 inch real color touch LCD for both local and remote (Remote is optional)	
ISI	Synchronization of injection delay and scanning delay	
Pressure limit	100 ~1200 Psi	100 ~600 Psi
Flow rate	0.1mL/s~50.0mL/s Step: 0.1mL/s for 0.1mL/s≤X≤10.0mL/s 1mL/s for 10mL/s<X≤50.0mL/s	0.1mL/s~25.0mL/s Step: 0.1mL/s for 0.1mL/s≤X≤10.0mL/s 1mL/s for 10mL/s<X≤25.0mL/s

SinoAngio-600/1200